

High dimensional morphology analysis reveals new insights in melanoma cell heterogeneity and enables label-free phenotyping

Evelyn Lattmann¹, Aizhan Tastanova¹, Andreja Jovic², **Manisha Ray²**, Tiffine Pham², Christian Corona², Jeanette Mei², Michael Phelan², Stephane C. Boutet², Ryan Carelli², Kevin B Jacobs², Julie Kim², Zhouyang Lian², Kiran Saini², Chassidy Johnson², Nianzhen Li², Mahyar Salek², Maddison Masaeli², Mitch Levesque¹
¹University Hospital Zurich, Schlieren, Switzerland. ²Deepcell Inc., Menlo Park, CA

KEY POINTS

- The Deepcell platform characterizes & sorts cells based on multi-dimensional morphology analysis without labels, eliminating the need for specific biomarkers.
- We demonstrate that multi-dimensional morphology can distinguish mesenchymal vs melanocytic cells in both cell lines and dissociated tumor biopsies
- An AI-based random forest classifier can predict cell phenotype based on morphology alone, verified by single-cell transcriptional information

INTRODUCTION

- Melanomas are the deadliest skin cancers, in part due to cellular plasticity and heterogeneity within the tumors.
- The Deepcell platform enables multi-dimensional morphology analysis and enrichment of unlabeled single cells using artificial intelligence (AI), advanced imaging, and microfluidics, enabling high resolution profiling of population heterogeneity.
- We imaged and analyzed 18 patient-derived melanoma cell lines representing both mesenchymal and melanocytic phenotypic states.
- High-dimensional morphological analysis showed distinct clusters for each phenotype, indicating distinct morphotype for each phenotype.
- We developed a random forest classifier to identify the top differential morphological features between the different cell lines, thereby providing a label-free means of phenotyping melanoma samples.
- The morphology analysis of the cell lines uncovered significant variability in pigmentation; a random forest classifier distinguished pigmented vs non-pigmented cells with >90% accuracy.
- Application of the classifier to images of dissociated tumor biopsies identified their phenotype, which was verified by scRNASeq, demonstrating the application of a label-free phenotyping using morphology alone in clinical samples.

METHODS




Figure 1. The Deepcell workflow. We analyzed each of the sample types using the Deepcell Workflow. **1)** Single cell suspension is **(2)** loaded onto a microfluidic chip. **(3)** Images of single cells are captured and analyzed in real-time by **(4)** deep learning and morphometric (computer vision) models to generate **(5)** multi-dimensional quantitative morphological profiles. User-defined cell clusters can be **(6a)** sorted for **(7)** downstream functional or molecular analysis. **(6b)** Morphology descriptions (embeddings) are extracted for data analysis/exploration and custom model training. The embeddings were used to profile the morphological heterogeneity of the cells, identify which morphology features differentiate each condition, and train a random forest classifier to predict each class of cells.

RESULTS




Figure 2. Multi-dimensional morphology distinguishes mesenchymal vs melanocytic cell lines. **(A)** 18 patient-derived cell lines from metastatic melanomas representing mesenchymal and melanocytic phenotypes were imaged and analyzed on the Deepcell platform. The phenotype of each cell line was based on mesenchymal or melanocytic gene expression signatures from bulk RNA-Seq data. **(B)** Images of cells indicate qualitative morphology differences detectable by eye, with melanocytic cells appearing smoother and mesenchymal cells are more granular. **(C)** Leiden clustering highlights the morphological diversity across the patient derived cell lines. **(D)** Multi-dimensional morphology UMAP showing distinct clusters for mesenchymal (orange) and melanocytic (green) cells, further illustrated by density plots. Inspection of the individual cell lines shows multiple morphotypes within each phenotype.

A random forest classifier can predict phenotype based on morphology alone




Figure 3. A random forest classifier can predict phenotype of melanoma cells with up to 77% accuracy. **(A)** We developed a random forest classifier to predict the phenotype of melanoma cells based on the images taken on the Deepcell platform. 8 cell lines with the most clear transcriptional score were used to train the classifier, which was then tested on an additional 12 cell lines. The classifier can predict the phenotype with up to 77% accuracy. **(B)** The classifier is concordant with transcriptional phenotype for 17/18 cell lines. The classifier provides single-cell phenotype information, noted in the percent of cells classified as each phenotype for each cell line. **(C)** Inspection of the images verified the predicted heterogeneity. **(D)** Local binary pattern (LBP) features, a metric of texture and granularity, are the top differential morphometric features distinguishing mesenchymal and melanocytic phenotypes. The distribution of the top feature is shown, suggesting mesenchymal cells (blue) have higher measures of granularity compared to melanocytic cells (orange). DL: deep learning.

RESULTS cont.

Morphometric analysis of morphology elucidates differences in pigmentation in phenotypic states




Figure 4. High dimensional morphology analysis reveals distinct pigmented subpopulation within melanoma cell lines. **(A)** Leiden clustering using morphometric features identified a population of pigmented cells (cluster 7). **(B)** Visual inspection and morphometric quantification of individual melanoma cell lines (MNT-1, 991104, SkMel-1) showed high, medium, and low pigmentation. **(C)** A random forest classifier using the morphometric features predict pigmented populations with up to 99% accuracy. **(D)** Top differential features distinguishing pigmented versus non-pigmented cells include pixel contrast and intensity values, as expressed by 25th/75th percentile of pixel grayscale values (Percentile_25/75), negative fraction, and mean pixel intensity.

scRNA-Seq verifies AI predictions of phenotype on melanoma DTCs




Figure 5. The phenotype classifier can predict the morphology of dissociated tumor biopsies. **(A)** The Melanoma Phenotype Classifier was run on dissociated tumor cell (DTC) samples, and AI predictions identified all DTCs as melanocytic. **(B)** We performed scRNA-Seq on the DTCs and the resulting UMAPs are shown. We analyzed the composition of the DTCs with sufficient malignant cell abundance for statistical analysis (labeled in blue). Results verified the AI predictions that the melanocytic phenotype was more prevalent in these samples.

CONCLUSIONS

- High dimensional morphology analysis shows multiple morphotypes for melanocytic and mesenchymal cells
- These morphological differences were used to develop a Melanoma Phenotype Classifier that can predict phenotype with up to 77% accuracy
- The morphometric features also identified a unique pigmented population with up to 99% accuracy
- The Melanoma Phenotype Classifier can be used to predict the phenotype of dissociated melanoma biopsies.
- These results suggest that high dimensional morphology can be used to characterize phenotype in a label-free manner, and provide new insights into tumor biology.

- Visit our other posters:
- Deepcell Platform [5381]
 - Deepcell AI & Data Science [5361]
 - Heterogeneous Tumor Subpopulation Enrichment [2392]
 - Malignant Effusion Tumor Cell Enrichment [LB170]

Learn more: www.deepcell.com

