# AACR Poster #1073 Al-assisted high-dimensional cytological single cell analysis distinguishes normal and malignant urine samples using cell morphology images with high accuracy at the single cell level

Thierry P.P. van den Bosch<sup>1\*</sup>, Geert J.L.H. van Leenders<sup>1\*</sup>, Siamala Sinnadurai<sup>1</sup>, Zhouyang Lian<sup>2</sup>, Ryan Carelli<sup>2</sup>, Matthew Nakaki<sup>2</sup>, Vivian Lu<sup>2</sup>, Janifer Cruz<sup>2</sup>, Kelly R. Miller<sup>2</sup>, Mahyar Salek<sup>2</sup>, Maddison M. Masaeli<sup>2</sup>, Peter J. van der Spek<sup>1</sup> <sup>1</sup> Erasmus MC Dept of Pathology & Clinical Bioinformatics Rotterdam, NL, <sup>2</sup> Deepcell Inc. Menlo Park CA, USA

### Introduction

Urothelial cancer (UCC) is the 4th most common malignancy in men. The current clinical problem in UCC management is the high interobserver variability of tissue grading. Moreover, low sensitivity of urine cytology for early detection of cancer, recurrences and disease progression remains a challenge.

Squamous-cell carcinoma (SCC) is a type of cancer that originates in squamous epithelial cells. Besides the epithelial lining of the skin, squamous cells also form the surface of the hollow organs that potentially are prone to viral infection. The organs most commonly affected by SCC are the skin, lungs, cervix and vulva, with rare occurrences (2-5%) also reported in the bladder.

High resolution images of single cells were captured using the Deepcell's REM-I platform. We trained a deep learning model on a cohort of urine samples consisting of 54 healthy controls, 19 atypical cases, and 14 bladder cancer cases. Applications of this study include improved early detection and quantitative grading of bladder cancer based on distinct morphological characteristics of cells. This deep learning model can potentially be re-used for other cancers derived from e.g. lung, cervix and vulva.

## **Study Goals**

- Label-free high-dimensional morphological profiling using Deepcell's REM-I to distinguish benign, atypical and malignant cells in urine cytology samples.
- Use Artificial Intelligence and Machine Learning for sorting and recognizing cells in urine, based on morphological differences.
- Count and collect the individual cells, belonging to distinct morphological clusters, to gain further insight in cell differentiation linked to bladder cancer.
- Sorted cells can be used for downstream analysis to study UCC/SCC relevant pathways, such as mutation detection (e.g. FGFR3), RNA expression profiling, immuno-histochemistry or proteomics.



Figure 1. The REM-I workflow. A single cell suspension is loaded onto a microfluidic chip where images of single cells are captured and analyzed in real-time by the Human Foundation Model (HFM). High-dimensional morphological features are visualized in a UMAP. User-defined cell clusters can be sorted for downstream functional or molecular analysis. Morphology embeddings are used to profile the morphological heterogeneity of the cells, identify which morphological features differentiate each condition, and train a random forest classifier.











